

# ReJoy™

FOR  
EXPANDABLE JOYSTICKS

ACTIVE CONTROLS  
REHAB SOLUTIONS


Position a  
JoyStick at  
Midline

Effortless Drop-In Connection • Robust MidLine Drive • Multiple Socket Locations

## ReJoy for Expandable Joysticks

The ReJoy facilitates the redeployment of most joysticks at midline. This is particularly important for the most compromised clients. In many instances, it is difficult or impossible for these clients to access an input device mounted on an armrest due to tone or ROM restrictions. Accessing an armrest mounted input device often requires the user to compromise their posture to drive a powerchair, which may result in permanent upper body deformations for active users.

The ReJoy mounting platform and lateral adjustable lower arm/hand gel support pads offers superior positioning of the input device which leads to increased function of the controls. The ReJoy actually works in conjunction with the seating system to maintain body position at midline, thereby maximizing client comfort while reducing fatigue.



Available Pre-Assembled • Interfaces with Most Expandable Electronics

Up to 3 Socket Positions • Gold Plated Connection Pins

Active Controls 597 Mantua Blvd. Sewell NJ 08080 USA

T 1.800.324.1527 E Sales@ActiveControls.com W ActiveControls.com

